

ANIMAL MEDICATION

SPECIAL ACID UNIT



Specially developed for use with concentrated propionic and formic acid products.



Includes PVDF body & special kalrez seals.



Suitable for pigs, poultry & other livestock.

D25RE09AO

0.1% to 0.9% dose rate10 to 2500 litres per hour flow rate0.3 to 6 bar working pressure

Temperature: 40°C maximum

Self Priming at Suction





Unit 4 Henson Park, Henson Way, Kettering, Northamptonshire, NN16 8PX tel. +44 (0)1536 461441 / e-mail: info@hingerose.co.uk www.hingerose.co.uk